














1) MOUNTING ARMS MAY BE REMOVED, AND POSITIONED ON EITHER SIDE OF THE CAGE PLATE AS NEEDED WITHOUT CHANGING THE ORIENTATION OF THE CAGE PLATE.
2) LARGER OPTICS MAY BE MOUNTED OUTSIDE OF THE CAGE SYSTEM WITH THE USE OF LONGER ER SERIES RODS (NOT INCLUDED).

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION		
	NAME	DATE
DRAWN	OA	14/JUL/21
APPROVAL	AK	15/JUL/21

COPYRIGHT © 2021 BY THORLABS

VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUNDOFF ERRORS

THORLABS
www.thorlabs.com

CYLINDRICAL LENS 30mm CAGE MOUNT

MATERIAL	REV
N/A	B

ITEM #	APPROX WEIGHT
CYCPA	0.06 kg